

The Biomechanics of Ankle-Foot Orthoses

Properly designed devices can positively affect gait.

Goals and Objectives

- 1) To provide an overview of the various mechanisms of action of ankle-foot orthoses on the lower extremity.
- 2) To compare the effects of solid versus articulated ankle-foot orthoses on the knee and midfoot joints.
- 3) To review research relevant to ankle braces and ankle-foot orthoses and their effects on balance and proprioception.
- 4) To provide evaluation, prescription and casting guidelines to optimize treatment outcomes with ankle-foot orthoses.

Welcome to Podiatry Management's CME Instructional program. Our journal has been approved as a sponsor of Continuing Medical Education by the Council on Podiatric Medical Education.

You may enroll: 1) on a per issue basis (at \$20.00 per topic) or 2) per year, for the special introductory rate of \$139 (you save \$61). You may submit the answer sheet, along with the other information requested, via mail, fax, or phone. In the near future, you may be able to submit via the Internet.

If you correctly answer seventy (70%) of the questions correctly, you will receive a certificate attesting to your earned credits. You will also receive a record of any incorrectly answered questions. If you score less than 70%, you can retake the test at no additional cost. A list of states currently honoring CPME approved credits is listed on pg. 182. Other than those entities currently accepting CPME-approved credit, Podiatry Management cannot guarantee that these CME credits will be acceptable by any state licensing agency, hospital, managed care organization or other entity. PM will, however, use its best efforts to ensure the widest acceptance of this program possible.

This instructional CME program is designed to supplement, NOT replace, existing CME seminars. The goal of this program is to advance the knowledge of practicing podiatrists. We will endeavor to publish high quality manuscripts by noted authors and researchers. If you have any questions or comments about this program, you can write or call us at: **Podiatry Management, P.O. Box 490, East Islip, NY 11730, (631) 563-1604 or e-mail us at bblock@podiatrym.com.**

Following this article, an answer sheet and full set of instructions are provided (p. 182).—**Editor**

By Douglas H. Richie Jr. DPM

Introduction

Prior to 1995, the dispensing of ankle-foot orthoses (AFOs) by podiatric physicians was a rare occur-

rence. Despite the fact that podiatric medical education provided stellar training in lower extremity biomechanics and taught the pathomechanics of virtually every neuromuscular disease, treatment op-

tions utilized by podiatric physicians in the '70's thru '90's were limited to foot orthotic therapy and shoe modifications.

In 1995, a modified version of

Continued on page 174

AFOs...

an ankle-foot orthosis, utilizing podiatric biomechanical principles (i.e., Richie Brace) was introduced to the profession. Since that time, the majority of podiatric physicians have expanded their treatment protocols to include AFO therapy into non-operative interventions for many lower extremity pathologies. These treatments have not been limited to podiatric devices, but now include almost every type of sophisticated brace, which were previously dispensed only by orthotists.

Bracing the lower extremity can be challenging and rewarding for the podiatric physician. However, evaluating patients with significant deformity and gait impairment requires the highest level of understanding of lower extremity biomechanics. This article is intended to provide an overview of the biomechanics of ankle-foot orthotic ther-

Figure 1: A three-point force system is applied by an ankle-foot orthosis. Control of joint moments by an orthosis requires a three-point fixation. This strategy applies three controlling forces to the limb: one placed over the joint center, with the other two forces applied proximal and distal to the joint and acting in opposite directions to each other.




Figure 2: Straps and pads applied to a podiatric ankle foot orthosis can control translational forces, i.e., internal/external rotation of the tibia

apy relevant to podiatric practice. Guidelines for evaluation, casting, and prescription of these devices will be provided to optimize outcomes with this powerful treatment intervention.

Definitions and Terminology

An orthosis is defined as: An externally applied device used to modify the structural or functional characteristics of the neuromusculo-skeletal system.¹ An orthosis has also been defined as: An apparatus used to support, align, prevent, or correct deformities or to improve the function of movable parts of the body.² Orthoses is the plural of orthosis.

The term "orthotic" has been traditionally used as an adjective, i.e., "orthotic device." Today, most dictionaries list both an adjective and a noun usage of the term "orthotic" and consider it synonymous with the term "orthosis."

Lower extremity orthoses function by affecting joint moments. A moment is also known as a force couple, which acts at a distance from an axis of rotation of a specific joint in the body. This force is measured in Newton-meters.

Joint moments can be produced externally and internally. External moments are produced by inertia and ground reaction forces. Internal moments are produced by active and passive structures within the body itself, such as muscles, ligaments and joint capsules. In the scientific literature, the term "moment" refers to internal moment, produced by anatomic structures which are actually resisting external moments which are

produced by physical forces. When internal structures fail, orthoses can modify external forces and moments to allow the body to function in a "normal" manner.

Internal moments are produced by active and passive structures within the body itself such as muscles, ligaments and joint capsules.

An external device used to support or improve function of the foot and ankle can take many physical forms. This orthotic can be as simple as a felt pad placed under the metatarsals or as sophisticated as a composite brace controlling foot and ankle motions.

Orthotics prescribed for lower extremity pathologies include foot orthoses (FOs), ankle-foot orthoses (AFOs), knee orthosis (KOs) and knee-ankle-foot orthoses (KAFOs). Podiatric physicians primarily prescribe or dispense FOs and AFOs.

Bowker has described four different ways in which an orthosis may modify the system of external forces and moment acting across a joint.³ The first three patterns of action are peculiar to ankle-foot orthoses, while the last pattern would also include the action of foot orthoses.

Restriction of Rotational Motion at a Joint

Restriction of rotational motion through control of joint moments is the most common reason for prescribing ankle foot orthoses in podiatric medicine. An AFO may limit the range of motion about any particular axis, or may limit the number of axes about which motion may occur. For example, an ankle foot orthosis can limit inversion/eversion motion of the ankle joint, while preserving dorsiflexion/plantarflexion.

Control of joint moments by an orthosis requires a three-point fixation. This strategy applies three controlling forces to the limb: one placed over the joint center, with the other two forces applied proximal and distal to the joint and act-

Continued on page 175

AFOs...

ing in opposite directions to each other (Figure 1). Applications of these forces proximal and distal to the subtalar and ankle joint are not possible with foot orthoses, but can be accomplished with ankle foot orthoses.

Ankle-foot orthoses can apply a moment to a joint not only to restrict motion, but also position a segment into a preferred alignment. For example, a fixed position solid AFO can maintain the foot at a pre-determined alignment at the ankle joint. For drop foot conditions, a solid AFO can hold the foot at a 90 degree angle to the leg. Clinicians must carefully evaluate available joint range of motion when prescribing AFO devices. If a patient has a fixed equinus deformity, fitting a neutral-positioned, 90 degree solid AFO will be impossible. Similarly, when there is a fixed varus alignment of the hindfoot, application of an AFO can only correct the foot position to the available range of eversion at the ankle and subtalar joint.

Restrict Translational Motion at a Joint

Since normal joint motion in the lower extremity is rotational, any translational motion could be considered abnormal. Translational motion occurs when there is loss of ligamentous integrity around the

joint, which can occur after acute trauma or as a result of a chronic degenerative process.

Braces which limit translational motion require four points of fixation. This requires a rigid brace which incorporates straps and pads applied tightly to the skin surface. This feature is incorporated into podiatric ankle foot orthoses which are designed to control excessive transverse plane rotation of the ankle joint which

Translational motion occurs when there is loss of ligamentous integrity around the joint, which can occur after acute trauma or as a result of a chronic degenerative process.

occurs in chronic ankle instability and adult-acquired flatfoot (Figure 2).

Control of Axial Forces Across a Joint

Axial loading of joints in the lower extremity can contribute to pain and disability when the articular cartilage is damaged, or the underlying bony structures become deformed. An ankle-foot orthosis can off-load a specific joint when it is properly designed to transfer load to other anatomic structures. An example of such a device is the patellar tendon-bearing orthosis.³

Control of Line of Action of Ground-Reaction Force

Ground-reaction forces pass through the foot from heel strike to toe off and create a line of action on every joint in the lower extremity. These ground-reaction forces will therefore create a moment about each of these joints. Orthoses can affect the alignment of ground reaction forces to change joint moments. This is the primary function of foot orthoses (Figure 4).




Figure 3: A patellar tendon weight-bearing brace

Foot orthoses cannot apply a three-point force system to the subtalar joint or the ankle joint, but can possibly do so at the midtarsal joint. Thus, kinematic studies of foot orthotic treatment effects have shown very minimal improvements of alignment of the rearfoot.⁴

Conversely, many studies have documented positive effects of foot orthoses to alter joint moments of the lower extremity.⁵⁻¹⁰ Munderman has shown that the most important feature of foot orthoses to positively affect lower joint moments is contouring the device to the shape of the foot.¹¹ Thus, ankle foot orthoses which contour the footplate can be expected to offer the positive effects on joint moments that are found with foot orthoses alone.

Studies of Ankle Foot Orthoses: Kinetics and Kinematics

Kinetic and Kinematic effects of ankle-foot orthoses have been extensively studied.¹²⁻¹⁶ However, most of this research has focused on the effects of ankle-foot orthoses on patients with neuromuscular conditions. Few studies have been published on the effects of ankle foot orthoses in healthy subjects, and

Continued on page 176




Figure 4: Foot orthoses can alter alignment of ground reaction forces to change moment to the subtalar joint

AFOs...

virtually no studies have been conducted on sport applications of these types of devices.

Kitaoka, et al. studied the kinetic and kinematic effects of three types of ankle foot orthoses in twenty healthy subjects walking over ground.¹⁷ In the frontal plane,




Figure 5: A solid AFO restricts ankle plantarflexion and causes knee flexion moment to bring the forefoot to the ground

all three orthoses (a solid AFO with footplate, solid AFO with heel portion only, and articulated AFO with footplate) significantly reduced maximal hindfoot inversion, but did not affect eversion. The solid ankle AFO design significantly reduced both plantarflexion and dorsiflexion of the ankle, while the articulated ankle AFO did not affect ankle sagittal plane motion compared to the unbraced condition. Midfoot motion was reduced with the articulated AFO, and increased with the solid AFO. Cadence was reduced with the solid AFO's. All three braces were associated with decreased aft and medial shear forces compared to the non-braced condition.

Radtka, et al. studied the kinetic and kinematic effects of solid and hinged (articulated) ankle foot orthoses on nineteen healthy subjects during stair locomotion.¹⁸ A unilateral hinged (articulated) ankle-foot

orthosis produced kinematic and kinetic effects which were similar to subjects wearing no orthosis. The unilateral solid ankle-foot orthosis produced more abnormal ankle joint angles, moments and powers and more proximal compensations at the knee, hip and pelvis than the hinged AFO during stair locomotion. Subjects wearing either orthosis

walked slower during stair locomotion compared to the non-braced condition.

Huang, et al. also compared solid and articulated ankle-foot orthoses in terms of restriction of motion at the ankle, hindfoot, and forefoot.¹⁹ The subjects all had arthritis of the ankle either post-traumatic or osteoarthritis. A solid ankle-foot orthosis restricted ankle motion better than an articulated ankle foot

orthosis, but also caused greater sagittal plane motion across the midfoot joints.

Hartsell and Spaulding measured passive resistive torque applied throughout inversion range of motion of the ankle in healthy subjects and those with chronically unstable ankles.²⁰ A hinged semi-rigid, non-custom ankle-foot orthosis demonstrated significant increased passive resistive inversion

torque forces and restricted overall inversion motion better than a lace up ankle brace.

What remains obscure is an un-

Excessive transverse plane rotation of the ankle joint occurs in chronic ankle instability and adult-acquired flatfoot.

derstanding of the optimal range and plane of motion controlled by an ankle orthosis to achieve a desired treatment effect. The studies just cited have conclusively shown that restriction of motion of any joint in the lower extremity will have negative effects in the neighboring joints, both proximal and distal.

Solid Versus Articulated Ankle Foot Orthoses: Clinical Considerations

The kinematic studies comparing solid versus articulated ankle-foot orthoses suggest that changes occur during ambulation which can have significant clinical consequences. Not only are forces or moments transmitted to neighboring joints, but also alterations of neuro-sensory feedback can adversely affect balance and proprioception.

When an orthosis eliminates motion at the ankle joint, significant adaptation must occur by the




Figure 6: A heel rocker can slow or reduce knee flexion moment when wearing a solid AFO

patient to continue ambulation. The primary effects of a fixed or "locked" ankle position during gait, achieved by wearing a solid AFO brace, are at the knee joint.

When the ankle is fixed at 90 degrees by an AFO, *Continued on page 177*

AFOs...

the contact phase of gait will be prolonged, as the ankle cannot plantarflex to bring the forefoot to the ground. This places a significant external knee flexion moment to force the tibia forward and bring the forefoot to the ground (Figure 5). The knee can be protected from this damaging force if the patient's shoe can be beveled in the heel (i.e., heel rocker) or if a pad is placed in the heel section of the shoe (Figure 6).²¹

A fixed-ankle AFO will also cause abnormal knee joint moment during mid-stance. As the leg passes forward over the foot, a lack of ankle joint dorsiflexion will prevent forward migration of the tibia which will then cause an external knee extension moment. This may cause pain or even hyperextension of the knee during late mid-stance. Hullin has shown that a rocker sole modifica-

cles, there is a tendency for delayed knee extension at the end of mid-stance. A solid AFO will extend the knee in patients with weak quadriceps and prevent the "drop-down" normally seen at the knee of these

patients during gait. A solid AFO can also provide knee flexion moment to prevent recurvatum deformity at the knee. By re-directing ground-reaction forces posterior to the center of rotation of the knee joint, a flexion moment, rather than extension moment, will occur. This is made possible by orienting the footplate alignment of




Figure 7: A solid AFO can direct ground reaction forces posterior to the knee axis of rotation to reduce a recurvatum deformity.

the AFO in slight dorsiflexion (Figure 7).

Balance and Proprioception

The effects of ankle braces on postural control have been extensively studied. Baier and Hopf studied 22 athletes with functional instability of the ankle joint compared to 22 healthy athletes.²³ A significant improvement of postural control, as evidenced by reduced mediolateral sway velocity, was found in the instability group when wearing both a rigid and semi-rigid stirrup ankle brace. However, other studies performed on both healthy subjects and on subjects with functional ankle instability have failed to show any improvements of postural control with the use of ankle braces.²⁴⁻²⁷

Similarly, studies of effects of ankle-foot orthoses on balance and proprioception do not provide consistent findings.²⁸ Yet, studies of treatment effects of these devices commonly attribute any positive findings to improvements in proprioception.^{29,30} Interestingly, studies of foot orthoses, compared to ankle foot orthoses, show

much more consistent improvements of balance and postural control.³¹

The reason that foot orthoses, rather than ankle foot orthoses, can improve postural control may relate to two mechanisms. First, the contoured shape of a foot orthosis has been speculated to improve the receptor field for sensory feedback in the proprioceptive loop between the foot and the spinal cord.³¹ Many traditional ankle braces and ankle foot orthoses do not have contoured footplates. Second, the limitation of ankle motion in restrictive AFO's and certain ankle braces can adversely affect other sensory pathways for proprioception.

When the ankle joint is locked, mechanoreceptors within the liga-

Continued on page 178




Figure 8: An AFO attached to a leg with tibial varum will invert the foot the same degree of varus deformity in the leg.

The most important feature of foot orthoses to positively affect lower joint moments is contouring the device to the shape of the foot.

tion of the shoe will decrease the knee extension moment caused by a fixed ankle AFO.²²

This same effect of a solid AFO on knee moments can be used in a positive way to improve gait in patients with muscle weakness. In patients with weak quadriceps mus-

AFOs...

ments cannot be stimulated to provide proprioceptive input.²⁸ Secondly, restriction of ankle motion can inhibit the stretch reflex of the leg muscles and tendons, which are considered the most important component of the somatosensory system for balance and postural control.^{28,31}

These factors become important when considering bracing an individual with compromised balance and postural control. Ironically, these patients may otherwise be the most likely candidates for AFO therapy in podiatric practice: elderly, diabetic, and neurologically impaired individuals.

In elderly patients, falls are the leading cause of accidental death.^{32,33} Studies have shown that patients with diabetes are fifteen times more likely to suffer catastrophic falls than their age-matched counterparts.³⁴ Whether due to sensory neuropathy of diabetes, or due to motor neuropathy from stroke, patients with neurologic impairment are at risk for falling.³⁵⁻³⁸

Lavery, et al. previously suggested that the efforts to off-load the neuropathic foot may have deleterious effects on postural control.³⁹ Studies of ankle-foot orthoses have shown that AFOs which restrict ankle motion will inhibit postural control and decrease velocity during gait.^{40,41} Cattaneo, et al. showed that AFOs would improve static balance in patients with multiple sclerosis, but would compromise dynamic balance during gait.⁴²

Evaluation and Prescription Guidelines

Patients should be carefully observed in gait before any joint range of motion and muscle testing are performed. Knee flexion stability should be evaluated. In cases of

recurvatum or excessive knee flexion during stance, a solid AFO should be prescribed.

A plantarflexed foot position at touchdown may be due to: fixed equinus of the ankle, spasticity of ankle plantarflexors or weakness of ankle dorsiflexors. A hinged AFO with dynamic assist ankle joints can be prescribed for only one of these conditions: weak ankle dorsiflexors with a flaccid drop foot deformity. In cases of spasticity, a hinged AFO with a plantarflexion stop is recommended.

During gait, and in static stance, balance and postural control should be evaluated. A modified Romberg test with eyes open will detect loss of postural control. This will commonly be found in patients with diabetes and with adult-acquired flatfoot deformity. Careful consideration should be

made about using solid AFO devices and gauntlet style braces which limit ankle joint motion in these patients. It is recommended that all patients with balance problems be prescribed a cane to use when wearing solid AFO devices. A cane has been demonstrated to normalize postural control in neuro-pathic patients.³⁷

Off-weight-bearing range of motion of the ankle and subtalar joints should be carefully evaluated

before casting and prescribing an AFO. If the ankle joint is not capable of dorsiflexing to neutral (90 degree alignment of foot to leg) with the knee extended, then a

neutral solid AFO should not be prescribed. The footplate of the solid AFO should be oriented to the maximum available dorsiflexion position when there is restricted ankle dorsiflexion. A heel lift should then be applied to the shoe or to the AFO to assure proper alignment of the leg and knee.

Restricted subtalar joint motion, particularly in the direction of eversion, may mandate modification of the AFO prescription. In cases of severe uncompensated rearfoot varus, commonly seen in Charcot Marie

Tooth disease, traditional AFO's may not be able to prevent lateral instability. Lateral wedging of the forefoot and rearfoot of the AFO footplate should be considered in these cases. Also, pedorthic modifications including a lateral shoe flare and/or lateral wedging of the midsole can improve stability in cases of reduced rearfoot eversion.

In terms of posting AFO braces, the practitioner should be aware that these additions will not only tilt the footplate, but will tilt the entire leg portion of the brace. Therefore, adding a 3 degree forefoot varus post to the footplate of the AFO will tilt the leg portion of the device 3 degrees inverted to the floor. If the patient does not have an existing tibial varum alignment of 3 degrees or more, posting the brace will cause an ill-fitting brace and potential rocking between the heel and forefoot.

At the same time, a patient with significant genu varum or tibial varum, exceeding 6 degrees, will require a modification of the AFO prescription. Unlike wearing foot orthoses where leg alignment will not change orthotic alignment, ankle-foot orthoses create a new situation not always familiar to the podiatric physician. The frontal plane alignment of the leg to the

Continued on page 179

When the ankle is fixed at 90 degrees by an AFO, the contact phase of gait will be prolonged, as the ankle cannot plantarflex to bring the forefoot to the ground.

The unilateral solid ankle-foot orthosis produces more abnormal ankle joint angles, moments and powers, and more proximal compensations at the knee, hip and pelvis than the hinged AFO during stair locomotion.

AFOs...

floor in midstance will determine the alignment of the entire AFO, including the footplate. In other words, the AFO, when attached to the leg, will align the footplate to the floor in the same alignment that the leg is oriented to the floor. Thus, when the tibia is aligned in 10 degrees of varum, attaching an AFO to this leg will orient the footplate 10 degrees inverted (Figure 8).

Significant tibial varum or genu varum will require that the footplate of the AFO be posted to the same degrees of limb varus to provide a stable interface between the AFO and the supportive surface. Alternatively, depending on the pathology, the footplate of the AFO can be hinged and re-oriented to rest flat on the supportive surface while the leg upright portion of the brace remains aligned in varus. This is recommended only when there is available range of motion in the ankle and subtalar joints to evert the foot away from the inverted tibia in order to bring the foot flat on the floor.

A partial weight-bearing cast, placing the foot flat on the floor, will capture any varus/valgus limb alignment which can be preserved in the fabrication of the AFO. However, in most cases, this technique will eliminate the ability of the practitioner to position the foot into subtalar neutral while the midtarsal joint is locked and stable.

The neutral suspension casting technique continues to be the preferred method of capturing an optimal model of the foot

for functional foot orthoses fabrication by podiatric physicians. Munderman has shown that the most important feature of foot orthotic devices to improve kinetics and kinematics of lower limb function is the contouring of the device to the shape of the foot.¹¹ Such contouring is lost when the impression casting

technique is performed on a weight-bearing foot. In a previous article published in *Podiatry Management* (Sept. 2007), the benefits of neutral suspension casting for AFO braces was demonstrated in comparison to partial and full weight-bearing casting techniques.

When neutral suspension casting techniques are utilized, leg alignment to the floor in static stance should be taken into consideration when balancing and fabricating the AFO brace.

Footwear should be considered and prescribed to match the AFO and patient needs. Solid AFO braces will translate significant abnormal moments to the knee both during contact phase and late mid-stance. A padded heel, a heel rocker, and a forefoot rocker are all shoe modification which can neutralize most of these abnormal knee moment forces. Solid AFOs and gauntlet AFOs will not fit into the posterior margin of the heel counter-insole

junction. This is because the solid plastic posterior heel of the brace will contact the top-line of the heel counter, which is always tilted slightly forward. This will push the brace forward in the shoe, and keep the foot plate from conforming to the posterior margin of the heel counter. Thus, solid AFOs and gauntlet braces

will mandate an increase of at least one full shoe size.

Finally, the shoe upper is an essential component of the three point force system which AFOs utilize to correct foot and leg alignment. The midfoot and forefoot shoe upper sections (vamp and toebox) should fit snugly and con-

form to the foot to prevent frontal, sagittal and transverse plane movement.

Summary

Ankle foot orthoses have significant advantages over foot orthoses to improve alignment and change joint moments in the lower extremity. However, application of AFOs can also have deleterious effects on a patient, which must be considered when prescribing these devices. Ankle-foot orthoses, depending on design, can have damaging effects on the knee and midfoot. These negative effects can be neutralized by appropriate shoe modifications.

The effects of bracing the lower extremity can also be positive or negative in terms of balance and proprioception. Many patients who require AFO therapy are already at risk for catastrophic falls. Practitioners should consider whether solid vs. articulated AFO's should be prescribed, taking into consideration the desired treatment effects versus the potential negative aspects of limiting motion around the ankle joint. ■

References

- 1 Dorlands Medical Dictionary, W.B. Saunders, 2000.
- 2 Bunch Wh, Keagy RD, Ritter A E, et al: Atlas of Orthotics: Biomechanical Principles and Application. 2nd Ed. C.V. Mosby, St. Louis, 1985; 151-159.
- 3 Bowker P, Condie DN, Bader DL, Pratt DJ. Biomechanical Basis of Orthotic Management. Butterworth Heinemann, Oxford, 1993; 27-38.
- 4 Nigg BM, Stergiou P, Cole G, Stefanyshyn D, Mundermann A, Humble N. Effect of shoe inserts on kinematics, center of pressure, and leg joint moments during running. 2003, Med Sci Sports Exerc 35(2): 314-9.
- 5 Stefanyshyn DJ, Stergiou P, Lun VMY, Meeuwisse WH, Nigg BM. Knee joint moments and patellofemoral pain syndrome in runners. Part I: A case control study; Part II: A prospective cohort

Continued on page 180

A padded heel, a heel rocker, and a forefoot rocker are all shoe modification which can neutralize most of these abnormal knee moment forces.

A solid AFO will extend the knee in patients with weak quadriceps and prevent the "drop-down" normally seen at the knee of these patients during gait.

AFOs...

study. *Int Soc Biomech: 4th Symposium on footwear biomechanics*, p 86.

⁶ Crenshaw SJ, Pollo FE, Calton EF. Effects of lateral-wedged insoles on kinetics at the knee. *Clin Orthop*, 2000, 375: 185-192.

⁷ Neptune RR, Wright IC, Van den Bogert AJ. The influence of orthotic devices and vastus medialis strength and timing on patellofemoral loads during running. *Clin Biomech* 2000, 15: 611-618.

⁸ Calton EF, Crenshaw SJ, Dill SM, Pollo. The effects of lateral wedged insole on kinetics and kinematics at the knee. 4th Ann. Gait Clin Movement Analysis Meeting, Dallas, Texas.

⁹ Nigg BM, Cole G, Stergiou P, Stefanyshyn D. The use of pressure measurements to determine the effect of shoe orthotics on knee joint moments. In: *Proceedings of the EMED Meeting*, July, Munich, Germany.

¹⁰ Williams DS, Davis IM, Baitch SP. Effect of inverted orthoses on lower-extremity mechanics in runners. *Med Sci Sport Exer*, 2003, 35(12): 2060-2068.

¹¹ Mundermann A, Nigg BM, Humble RN, Stefanyshyn DJ. Foot orthotics affect lower extremity kinematics and kinetics during running. *Clin Biomech* 2003, 18(3):254-262.

¹² Romkes J, Brunner R. Comparison of a dynamic and a hinged ankle-foot orthosis by gait analysis in patients with hemiplegic cerebral palsy. *Gait Posture* 2002; 15:18-24.

¹³ Franceschini M, Massucci M, Ferrari L, Agosti M, Paroli C. Effects of an ankle-foot orthosis on spatiotemporal parameters and energy cost of hemiparetic gait. *Clin Rehabil* 2003; 17:368-372.

¹⁴ Wang RY, Yen L, Lee CC, Lin PY, Wang MF, Yang YR. Effects of an ankle-foot orthosis on walking ability in chronic stroke patients: a randomized controlled trial. *Clin Rehabil* 2004; 18:550-557.

¹⁵ Geboers JF, Drost MR, Spaans F, Kuipers H, Seelen HA. Immediate and long-term effects of ankle-foot orthosis on muscle activity during walking: a randomized study of patients with unilateral foot drop. *Arch Phys Med Rehabil* 2002; 83:240-245.

¹⁶ de Wit DC, Buurke JH, Nijlant JM, Izerman MJ, Hermens HJ. The effect of an ankle-foot-orthosis on walking ability in chronic stroke patients: a randomized controlled trial. *Clin Rehabil* 2004; 18:550-557.

¹⁷ Kitaoka HB, Crevoisier XM, Harbst K, Hansen D, Kotajarvi B, Kaufman K. The effect of custom-made braces for the ankle and hindfoot on ankle and

foot kinematics and ground reaction forces. *Arch Phys Med Rehabil* 2006; 87:130-135.

¹⁸ Radtka SA, Oliveira GB, Lindstrom KE, Borders MD. The kinematic and kinetic effects of solid, hinged, and no ankle-foot orthoses on stair locomotion in healthy adults. *Gait Posture* 2005; (2):211-8.

¹⁹ Huang YC, Harbst K, Kotajarvi B, Hansen D et al. Effects of ankle-foot orthoses on ankle and foot kinematics in patients with ankle osteoarthritis. *Arch Phys Med Rehabil* 87: 710-716, 2006.

²⁰ Hartsell HD, Spaulding S J. Effectiveness of external orthotic support on passive soft tissue resistance of the chronically unstable ankle. *Foot Ankle Int* 1997; 18:144-150.

²¹ Wiest DR, et al.: The influence of heel design on a rigid ankle foot orthosis. *Orthot Prosthetics* 33:3, 1979

²² Hullin MG, Robb JE. Biomechanical effects of rockers on walking in a plaster

²³ Baier M, Hopf T. Ankle orthoses effect on single-limb standing balance in athletes with functional ankle instability. *Arch Phys Med Rehabil* 1998; 79:939-944

²⁴ Kinzey SJ, Ingersoll CD, Knight KL. The effects of selected ankle appliances on postural control. *Jour Athl Train* 1997; 32:300-303.

²⁵ Palmieri RM, Ingersoll CD, Cordova ML, Kinzey SJ. The spectral qualities of postural control are unaffected by 4 days of ankle-brace application. *J Athl Train* 2002; 37:269-274.

²⁶ Barkoukis V, Sykaras E, Costa F, Tsoarbatzoudis H. Effectiveness of taping and bracing in balance. *Percept Mot Skills* 2002; 94:566-574.

²⁷ Wikstrom EA, Arrigenna MA, Tillman MD, Paul AB. Dynamic postural stability in subjects with braced, functionally unstable ankles. *J Athl Train* 2006; 41:245-250.

²⁸ Hijmans J M, Geertzen JAB, Dijkstra PU, Postema K. A systematic review of the effects of shoes and other ankle or foot appliances on balance in older people and people with peripheral nervous system disorders. *Gait and Posture* 25: 316-323, 2007.

²⁹ Sittler M, Ryan J, Wheeler B, MCBride J, Arciero R, Anderson J, Horodyski M. The efficacy of a semirigid ankle stabilizer to reduce acute ankle injuries in basketball: a randomized clinical study at west Point. *Am J Sports Med*. 1994; 22:454-461.

³⁰ Surve I, Schweltnus MP, Noakes T, Lombard C. A fivefold reduction in the incidence of recurrent ankle sprains in soccer players using the sport-stirrup orthosis. *Am J Sports Med*. 1994; 12:601-606.

³¹ Richie DH. Effects of foot or-

thoses to treat chronic ankle instability. *Jour Am Pod Med Assoc* 97(1): 19-30, 2007.

³² Stel VS, Smit JH, Pluijms SMF, Lips P. Balance and mobility performance as treatable risk factors for recurrent falling in older persons. *J Clin Epidemiol* 56: 659-669, 2003.

³³ Clark S, Rose DJ, Fujimoto K. Generalizability of the limits of stability test in the evaluation of dynamic balance among older adults. *Archives of Physical Med Rehabil* 78: 587-591, 1997.

³⁴ Cavanagh PR, Derr JA, Ulbrecht JS et al. Problems with gait and posture in neuropathic patients with insulin-dependent diabetes mellitus. *Diabetic Med* 9: 469-474, 1992.

³⁵ Richardson JK, Ashton-Miller JA. Peripheral neuropathy. An overlooked cause of falls in the elderly. *Postgraduate Medicine* 99: 161-172, 1996.

³⁶ Simoneau GG, Ulbrecht JS, Derr JA, Becker MB, Cavanagh PR. Postural instability in patients with diabetic sensory neuropathy. *Diabetes Care* 17: 1411-1421, 1994.

³⁷ Richardson JK, Hurvitz EA. Peripheral neuropathy: a true risk factor for falls. *J Gerontol A Biol Sci Med Sci* 50: M211-15, 1995.

³⁸ Simmons RW, Richardson C, Pozos R. Postural stability of diabetic patients with and without cutaneous sensory deficit in the foot. *Diabetes Res Clin Pract* 36: 153-160, 1997.

³⁹ Lavery LA, Fleishl JG, Laughlin TJ, Vela SA, Lavery DC, Armstrong DG. Is postural stability exacerbated by off-loading devices in high risk diabetics with foot ulcers? *Ostomy Wound Manage* 44: 26-34, 1998.

⁴⁰ Rao N, Aruin A. The effect of ankle-foot-orthoses on balance impairment: single case study. *Jour Prost Orth* 1999; 11:15-21.

⁴¹ Nahornizk MT et al. Kinematic compensations as children reciprocally ascend and descend stairs with unilateral and bilateral solid ankle-foot orthoses. *Gait Posture* 9: 199-206, 1999.

⁴² Cattaneo D, Marazzini F, Crippa A, Cardini R. Do static or dynamic AFOs improve balance? *Clinical Rehab* 2002.

Dr. Richie is Adjunct Associate Clinical Professor, Department of Applied Biomechanics, California School of Podiatric Medicine at Samuel Merritt College. He can be reached at drichejr@aol.com



See answer sheet on page 183.

- | | | |
|---|---|---|
| <p>1) An orthotic is:
 A) A metatarsal pad
 B) A custom foot support
 C) A pre-fabricated foot support
 D) All of the above</p> <p>2) Internal joint moments are produced by:
 A) Ground-reaction forces
 B) Inertia
 C) Cartilage
 D) Tendons and ligaments</p> <p>3) AFOs can limit motion better than foot orthoses because they:
 A) Have stronger plastic
 B) Have special hinges
 C) Apply three-point force systems
 D) Have special straps</p> <p>4) AFOs can decrease the pain of osteoarthritis of the knee by
 A) Off-loading compressive forces
 B) Cushioning the bottom of the foot
 C) Locking the ankle
 D) Assisting ankle joint motion</p> <p>5) Abnormal translational joint motion occurs when:
 A) Muscles become weak
 B) Ligaments are disrupted around the joint
 C) Proprioception is lost
 D) Internal joint moments increase</p> | <p>6) Abnormal transverse joint motion is seen in:
 A) Cerebral palsy
 B) Muscular Dystrophy
 C) Charcot Marie Tooth Disease
 D) Adult-Acquired Flatfoot</p> <p>7) Foot orthoses and AFOs can both:
 A) Alter the line of ground-reaction forces
 B) Improve skeletal alignment
 C) Reposition the ankle
 D) Lock the knee</p> <p>8) The most important feature of foot orthoses to affect joint moments is:
 A) Posting
 B) Molding or contouring the footplate
 C) Topcover
 D) Heel cup</p> <p>9) Studies have shown that mid-foot motion increases with:
 A) Dynamic assist AFOs
 B) Hinged AFOs
 C) Solid AFOs
 D) Patellar tendon bearing AFOs</p> <p>10) Proximal compensation at the knee, hip and pelvis is seen with:
 A) Solid AFOs
 B) Dynamic assist AFOs
 C) Hinged AFOs
 D) Articulated AFOs</p> | <p>11) Which phase of gait is prolonged when wearing a solid AFO?
 A) Midstance
 B) Heel rise
 C) Terminal stance
 D) Contact phase</p> <p>12) Solid AFOs cause a significant external moment causing:
 A) Knee extension at contact
 B) Knee flexion at contact
 C) Toe grasping
 D) Calf spasm</p> <p>13) During midstance, solid AFOs can cause:
 A) Pelvic tilt
 B) Knee flexion
 C) Knee hyperextension
 D) Ankle plantarflexion</p> <p>14) The negative effects of solid AFOs at the knee can be prevented by:
 A) Rocker shoe soles
 B) Special straps
 C) Padded topcovers
 D) Stronger plastic</p> <p>15) Patients with weak quadriceps muscles should be prescribed a:
 A) Hinged AFO
 B) Solid AFO
 C) Dynamic Assist AFO
 D) Custom foot orthotic</p> |
|---|---|---|

Continued on page 182

PM's CPME Program

Welcome to the innovative Continuing Education Program brought to you by *Podiatry Management Magazine*. Our journal has been approved as a sponsor of Continuing Medical Education by the Council on Podiatric Medical Education.

Now it's even easier and more convenient to enroll in PM's CE program!

You can now enroll at any time during the year and submit eligible exams at any time during your enrollment period.

PM enrollees are entitled to submit ten exams published during their consecutive, twelve-month enrollment period. Your enrollment period begins with the month payment is received. For example, if your payment is received on September 1, 2006, your enrollment is valid through August 31, 2007.

If you're not enrolled, you may also submit any exam(s) published in PM magazine within the past twelve months. **CME articles and examination questions from past issues of *Podiatry Management* can be found on the Internet at <http://www.podiatrym.com/cme>.** Each lesson is approved for 1.5 hours continuing education contact hours. Please read the testing, grading and payment instructions to decide which method of participation is best for you.

Please call (631) 563-1604 if you have any questions. A personal operator will be happy to assist you.

Each of the 10 lessons will count as 1.5 credits; thus a maximum of 15 CME credits may be earned during any 12-month period. You may select any 10 in a 24-month period.

The Podiatry Management Magazine CME program is approved by the Council on Podiatric Education in all states where credits in instructional media are accepted. This article is approved for 1.5 Continuing Education Contact Hours (or 0.15 CEU's) for each examination successfully completed.

16) If the footplate of a solid AFO is positioned in dorsiflexion, the effect on the knee in static stance is:

- A) Flexion
- B) Extension
- C) Varus
- D) Valgus

17) Drop foot can be controlled by a:

- A) Foot orthosis
- B) Rearfoot post
- C) Standard hinged AFO
- D) Solid AFO

18) A patient with midfoot arthritis can have more motion and pain when wearing a:

- A) Hinged AFO
- B) Solid AFO
- C) Dynamic Assist AFO
- D) Foot Orthosis

19) Balance and Proprioception can be compromised by AFOs due to:

- A) Restriction of joint range of motion
- B) Restrict input of muscle spindles
- C) Restrict input of joint proprioceptors
- D) All of the above

20) Tibial varum deformity of 10 degrees will orient the footplate of a neutral AFO:

- A) Flat on the floor
- B) 5 degrees inverted to the floor
- C) 10 degrees inverted to the floor
- D) 20 degrees inverted to the floor

See answer sheet on page 183.

**Home Study CME credits now
accepted in Pennsylvania**

Enrollment/Testing Information and Answer Sheet

Note: If you are mailing your answer sheet, you must complete all info. on the front and back of this page and mail with your credit card information to: **Podiatry Management, P.O. Box 490, East Islip, NY 11730.**

TESTING, GRADING AND PAYMENT INSTRUCTIONS

(1) Each participant achieving a passing grade of 70% or higher on any examination will receive an official computer form stating the number of CE credits earned. This form should be safeguarded and may be used as documentation of credits earned.

(2) Participants receiving a failing grade on any exam will be notified and permitted to take one re-examination at no extra cost.

(3) All answers should be recorded on the answer form below. For each question, decide which choice is the best answer, and circle the letter representing your choice.

(4) Complete all other information on the front and back of this page.

(5) Choose one out of the 3 options for testgrading: mail-in, fax, or phone. To select the type of service that best suits your needs, please read the following section, "Test Grading Options".

TEST GRADING OPTIONS

Mail-In Grading

To receive your CME certificate, complete all information and mail with your credit card information to:

**Podiatry Management
P.O. Box 490, East Islip, NY 11730**

There is **no charge** for the mail-in service if you have already enrolled in the annual exam CPME program, and we receive this

exam during your current enrollment period. If you are not enrolled, please send \$20.00 per exam, or \$139 to cover all 10 exams (thus saving \$61* over the cost of 10 individual exam fees).

Facsimile Grading

To receive your CPME certificate, complete all information and fax 24 hours a day to 1-631-563-1907. Your CPME certificate will be dated and mailed within 48 hours. This service is available for \$2.50 per exam if you are currently enrolled in the annual 10-exam CPME program (and this exam falls within your enrollment period), and can be charged to your Visa, MasterCard, or American Express.

If you are *not* enrolled in the annual 10-exam CPME program, the fee is \$20 per exam.

Phone-In Grading

You may also complete your exam by using the toll-free service. Call 1-800-232-4422 from 10 a.m. to 5 p.m. EST, Monday through Friday. Your CPME certificate will be dated the same day you call and mailed within 48 hours. There is a \$2.50 charge for this service if you are currently enrolled in the annual 10-exam CPME program (and this exam falls within your enrollment period), and this fee can be charged to your Visa, Mastercard, American Express, or Discover. If you are not currently enrolled, the fee is \$20 per exam. When you call, please have ready:

1. Program number (Month and Year)
2. The answers to the test
3. Your social security number
4. Credit card information

In the event you require additional CPME information, please contact PMS, Inc., at **1-631-563-1604**.

ENROLLMENT FORM & ANSWER SHEET

Please print clearly...Certificate will be issued from information below.

Name _____ Soc. Sec. # _____
Please Print: FIRST MI LAST

Address _____

City _____ State _____ Zip _____

Charge to: Visa MasterCard American Express

Card # _____ Exp. Date _____

Note: Credit card is the only method of payment. Checks are no longer accepted.

Signature _____ Soc. Sec.# _____ Daytime Phone _____

State License(s) _____ Is this a new address? Yes _____ No _____

Check one: I am currently enrolled. (If faxing or phoning in your answer form please note that \$2.50 will be charged to your credit card.)

I am not enrolled. Enclosed is my credit card information. Please charge my credit card \$20.00 for each exam submitted. (plus \$2.50 for each exam if submitting by fax or phone).

I am not enrolled and I wish to enroll for 10 courses at \$139.00 (thus saving me \$61 over the cost of 10 individual exam fees). I understand there will be an additional fee of \$2.50 for any exam I wish to submit via fax or phone.



**EXAM #7/09
The Biomechanics
of Ankle-Foot Orthoses
(Richie)**

Circle:

- | | |
|-------------|-------------|
| 1. A B C D | 11. A B C D |
| 2. A B C D | 12. A B C D |
| 3. A B C D | 13. A B C D |
| 4. A B C D | 14. A B C D |
| 5. A B C D | 15. A B C D |
| 6. A B C D | 16. A B C D |
| 7. A B C D | 17. A B C D |
| 8. A B C D | 18. A B C D |
| 9. A B C D | 19. A B C D |
| 10. A B C D | 20. A B C D |

LESSON EVALUATION

Please indicate the date you completed this exam

How much time did it take you to complete the lesson?

_____ hours _____ minutes

How well did this lesson achieve its educational objectives?

_____ Very well _____ Well

_____ Somewhat _____ Not at all

What overall grade would you assign this lesson?

A B C D

Degree _____

Additional comments and suggestions for future exams:

