

Reliable Drills from European Footcare Supply

European Footcare has been a partner to DPMs delivering routine footcare for more than 15 years. We have learned that podiatrists value their drills for powerful performance at the bur tip, excellent



Micro-Air

control, mitigation of dust, and reliability. That is why many of our users are still using drills that are more than 10 years old and why they contact European Footcare Supply when they require additional drills. New customers are surprised by the torque delivered to the bur tip that allows them to be more efficient in debridement of nail and callus. Dust is mitigated with an excellent vacuum that is surprisingly quiet. The vacuum nozzle evacu-

ates dust right at the bur tip where dust is created.

The Micro-Air vacuum drill has a serious medical appearance that patients recognize and are comfortable with. The vacuum hose is flexible and long enough for most clinics and has lifelong warranty. The handpieces are all metallic (no plastic here) and have a reliable locking chuck for safety and for easy exchange of the bur. All standard podiatry burs will lock in reliably in this handpiece.

The Ortho-Spray drill is a treatment modality that cools bur and patient tissue as you work. Dust binds with the water/alcohol spray and does not become airborne. Patients tolerate this treatment very well and the user



Orthospray

can work faster since there is no waiting for the heat in patient tissue to dissipate. Especially in small treatment areas such as with ingrown nail, the water spray drill is very efficient. If you have a laser for nail treatment, then thinning the nail plate prior to the laser application is

**After 2 weeks you will notice
that your treatment times are shorter
and you can see more patients.**

essential and patients appreciate the highly professional water spray treatment they received.

If this is the first time you're considering purchasing a dependable drill for your practice then insist on a trial period when you order your first drill. It is very important that you become comfortable with the ergonomic handpiece and the choice of burs available to you. You can now handle the carbide burs with confidence because of the better control this handpiece gives you. After 2 weeks you will notice that your treatment times are shorter and you can see more patients. That is why all high volume clinics use powerful drills.

Our drills come complete with footswitch, a cleaning tool and bur ready to go out of the box. The drills have 3 years warranty (that includes the handpiece). If you would like to discuss your needs for a drill or burs or German Stainless Steel hand instruments, *contact us at 1-800-881-2243, visit www.a-foot.com, or click here.* We are here for you.